

Key: >>> ... marks correct (expected) result

> ... marks conditionally correct (acceptable) result

± ... marks the result not evaluated

N ... number of the results

Nrel ... relative number of the results

Test

Sample	Group	N	Result	N _{rel}
--------	-------	---	--------	------------------

Set of samples 1

(575) Anti-thyreoglobulin

A1	(1) Abbott	>>>	49	Negative	28 %
	(12) Beckman Coulter	>>>	23	Negative	13 %
	(15) bioMérieux	>>>	1	Negative	0,58 %
	(25) TOSOH	>>>	4	Negative	2,3 %
	(29) Siemens (DPC)	>>>	19	Negative	11 %
	(36) Immunotech	>>>	3	Negative	1,7 %
	(60) Roche	>>>	54	Negative	31 %
	(135) OrgenTec	>>>	1	Negative	0,58 %
	(164) DiaSorin	>>>	2	Negative	1,2 %
	(179) Siemens (Bayer)	>>>	15	Negative	8,7 %
	(998) home-made	>>>	1	Negative	0,58 %
	(999) another manufacturer	>>>	1	Negative	0,58 %
----- Total per sample: n = 173		-----	Success: 100 %	-----	
B1	(1) Abbott	>>>	49	Positive	28 %
	(12) Beckman Coulter	>>>	23	Positive	13 %
	(15) bioMérieux	>>>	1	Positive	0,58 %
	(25) TOSOH	>>>	4	Positive	2,3 %
	(29) Siemens (DPC)	>>>	19	Positive	11 %
	(36) Immunotech	>>>	3	Positive	1,7 %
	(60) Roche		3	Negative	1,7 %
			1	Inconclusive result (borderline)	0,58 %
		>>>	50	Positive	29 %
	(135) OrgenTec	>>>	1	Positive	0,58 %
	(164) DiaSorin	>>>	2	Positive	1,2 %
	(179) Siemens (Bayer)	>>>	15	Positive	8,7 %
	(998) home-made		1	Negative	0,58 %
	(999) another manufacturer	>>>	1	Positive	0,58 %
----- Total per sample: n = 173		-----	Success: 97 %	-----	
Total success: 97 %					

(576) Anti-microsomes (TPO)

A1	(1) Abbott	>>>	51	Negative	28 %
	(12) Beckman Coulter	>>>	26	Negative	14 %
	(15) bioMérieux	>>>	1	Negative	0,56 %
	(25) TOSOH	>>>	5	Negative	2,8 %
	(29) Siemens (DPC)	>>>	20	Negative	11 %
	(36) Immunotech	>>>	2	Negative	1,1 %
	(60) Roche	>>>	52	Negative	29 %
			3	Positive	1,7 %
	(135) OrgenTec	>>>	1	Negative	0,56 %
	(164) DiaSorin	>>>	1	Negative	0,56 %
	(179) Siemens (Bayer)	>>>	16	Negative	8,9 %
	(998) home-made	>>>	1	Negative	0,56 %
	(999) another manufacturer	>>>	1	Negative	0,56 %
----- Total per sample: n = 180		-----	Success: 98 %	-----	
B1	(1) Abbott	>>>	51	Positive	28 %
	(12) Beckman Coulter	>>>	26	Positive	14 %
	(15) bioMérieux	>>>	1	Positive	0,56 %
	(25) TOSOH	>>>	5	Positive	2,8 %
	(29) Siemens (DPC)	>>>	20	Positive	11 %
	(36) Immunotech	>>>	2	Positive	1,1 %
	(60) Roche		1	Negative	0,56 %
		>>>	54	Positive	30 %
	(135) OrgenTec	>>>	1	Positive	0,56 %
	(164) DiaSorin	>>>	1	Positive	0,56 %
	(179) Siemens (Bayer)	>>>	16	Positive	8,9 %
	(998) home-made	>>>	1	Positive	0,56 %
	(999) another manufacturer	>>>	1	Positive	0,56 %
----- Total per sample: n = 180		-----	Success: 99 %	-----	
Total success: 98 %					

Test

Sample	Group	N	Result	N _{rel}
--------	-------	---	--------	------------------

Set of samples 2**(590) Anti-TSH receptor**

A2	(24) B.R.A.H.M.S.	>>>	6	Negative	6,5 %	
			1	Inconclusive result (borderline)	1,1 %	
	(29) Siemens (DPC)	>>>	9	Negative	9,7 %	
	(36) Immunotech	>>>	2	Negative	2,2 %	
	(60) Roche	>>>	71	Negative	76 %	
			1	Inconclusive result (borderline)	1,1 %	
	(87) Euroimmun	>>>	1	Negative	1,1 %	
	(146) MEDIPAN Diagnostica	>>>	1	Negative	1,1 %	
	(149) Siemens (Dade)	>>>	1	Negative	1,1 %	
	Total per sample: n = 93		Success: 98 %			
	B2	(24) B.R.A.H.M.S.	>>>	7	Positive	7,5 %
				9	Positive	9,7 %
		(29) Siemens (DPC)	>>>	2	Positive	2,2 %
(36) Immunotech		>>>	72	Positive	77 %	
(60) Roche		>>>	1	Positive	1,1 %	
(87) Euroimmun		>>>	1	Positive	1,1 %	
(146) MEDIPAN Diagnostica		>>>	1	Positive	1,1 %	
(149) Siemens (Dade)	>>>	1	Positive	1,1 %		
Total per sample: n = 93		Success: 100 %				
Total success: 98 %						