

Key:	>>> ... marks correct (expected) result	N ... number of the results
	> ... marks conditionally correct (acceptable) result	Nrel ... relative number of the results
	± ... marks the result not evaluated	

Test	Sample	Group	N	Result	N _{rel}
(472) anti-GAD					
A	(36)	Immunotech	>>>	1 Positive	3,0 %
	(49)	BioVendor	>>>	3 Positive	9,1 %
	(52)	Dynex	>>>	1 Positive	3,0 %
	(66)	DRG	>>>	5 Positive	15 %
	(87)	Euroimmun	>>>	15 Positive	45 %
	(146)	MEDIPAN Diagnostica	>>>	2 Positive	6,1 %
	(999)	another manufacturer	>>>	1 Inconclusive result (borderline)	3,0 %
			>>>	5 Positive	15 %
-----		Total per sample: n = 33	-----	Success: 97 %	-----
B	(36)	Immunotech	>>>	1 Positive	3,0 %
	(49)	BioVendor	>>>	3 Positive	9,1 %
	(52)	Dynex	>>>	1 Positive	3,0 %
	(66)	DRG	>>>	5 Positive	15 %
	(87)	Euroimmun	>>>	15 Positive	45 %
	(146)	MEDIPAN Diagnostica	>>>	2 Positive	6,1 %
	(999)	another manufacturer	>>>	6 Positive	18 %
-----		Total per sample: n = 33	-----	Success: 100 %	-----
Total success: 97 %					

(474) anti-IA-2					
A	(36)	Immunotech	>>>	1 Positive	3,2 %
	(49)	BioVendor	>>>	1 Negative	3,2 %
			>>>	2 Positive	6,5 %
	(52)	Dynex	>>>	1 Positive	3,2 %
	(66)	DRG	>>>	3 Positive	9,7 %
	(87)	Euroimmun	>>>	1 Negative	3,2 %
			>>>	12 Positive	39 %
	(146)	MEDIPAN Diagnostica	>>>	2 Positive	6,5 %
	(999)	another manufacturer	>>>	8 Positive	26 %
-----		Total per sample: n = 31	-----	Success: 94 %	-----
B	(36)	Immunotech	>>>	1 Negative	3,2 %
	(49)	BioVendor	>>>	2 Negative	6,5 %
			>>>	1 Positive	3,2 %
	(52)	Dynex	>>>	1 Negative	3,2 %
	(66)	DRG	>>>	3 Negative	9,7 %
	(87)	Euroimmun	>>>	12 Negative	39 %
			>>>	1 Positive	3,2 %
	(146)	MEDIPAN Diagnostica	>>>	2 Negative	6,5 %
	(999)	another manufacturer	>>>	8 Negative	26 %
-----		Total per sample: n = 31	-----	Success: 94 %	-----
Total success: 94 %					