

Key:	>>> ... marks correct (expected) result	N ... number of the results
	> ... marks conditionally correct (acceptable) result	Nrel ... relative number of the results
	± ... marks the result not evaluated	

Test

Sample	Group		N	Result	N _{rel}
(268) Bilirubin					
A	(5) Dialab	>>>	1	Negative	1,7 %
	(12) Beckman Coulter	>>>	2	Negative	3,3 %
	(18) Cormay	>>>	4	Negative	6,7 %
	(46) Erba Lachema	>>>	10	Negative	17 %
	(60) Roche	>>>	2	Negative	3,3 %
	(63) Sysmex	>>>	2	Negative	3,3 %
	(109) Axis-Shield	>>>	2	Negative	3,3 %
	(150) Arkray	>>>	2	Negative	3,3 %
	(161) YD Diagnostics	>>>	1	Negative	1,7 %
	(162) Siemens (Atellica)	>>>	3	Negative	5,0 %
	(179) Siemens	>>>	4	Negative	6,7 %
	(197) Analyticon	>>>	1	Negative	1,7 %
	(201) Iris Diagnostics	>>>	1	Negative	1,7 %
	(202) Dirui	>>>	19	Negative	32 %
	(207) 77 Elektronika Kft.	>>>	5	Negative	8,3 %
	(999) another manufacturer	>>>	1	Negative	1,7 %
Total per sample: n = 60		Success: 100 %			
B	(5) Dialab	>>>	1	3+, 4+	1,7 %
	(12) Beckman Coulter	>	2	1+, 2+	3,3 %
	(18) Cormay	>>>	4	3+, 4+	6,7 %
	(46) Erba Lachema	>>>	10	3+, 4+	17 %
	(60) Roche	>>>	2	3+, 4+	3,3 %
	(63) Sysmex	>>>	2	3+, 4+	3,3 %
	(109) Axis-Shield	>>>	2	3+, 4+	3,3 %
	(150) Arkray	>>>	2	3+, 4+	3,3 %
	(161) YD Diagnostics	>>>	1	3+, 4+	1,7 %
	(162) Siemens (Atellica)	>>>	3	3+, 4+	5,0 %
	(179) Siemens	>>>	4	3+, 4+	6,7 %
	(197) Analyticon	>>>	1	3+, 4+	1,7 %
	(201) Iris Diagnostics	>>>	1	3+, 4+	1,7 %
	(202) Dirui	>	1	1+, 2+	1,7 %
		>>>	18	3+, 4+	30 %
	(207) 77 Elektronika Kft.	>>>	5	3+, 4+	8,3 %
	(999) another manufacturer	>>>	1	3+, 4+	1,7 %
Total per sample: n = 60		Success: 100 %			
Total success: 100 %					

(264) Total protein

A	(5) Dialab	>>>	1	Negative, < 0,3 g/L	1,6 %
	(12) Beckman Coulter	>>>	2	Negative, < 0,3 g/L	3,3 %
	(18) Cormay	>>>	4	Negative, < 0,3 g/L	6,6 %
	(46) Erba Lachema	>>>	11	Negative, < 0,3 g/L	18 %
	(60) Roche	>>>	2	Negative, < 0,3 g/L	3,3 %
	(63) Sysmex	>>>	2	Negative, < 0,3 g/L	3,3 %
	(109) Axis-Shield	>>>	2	Negative, < 0,3 g/L	3,3 %
	(150) Arkray	>>>	2	Negative, < 0,3 g/L	3,3 %
	(161) YD Diagnostics	>>>	1	Negative, < 0,3 g/L	1,6 %
	(162) Siemens (Atellica)	>>>	3	Negative, < 0,3 g/L	4,9 %
	(179) Siemens	>>>	4	Negative, < 0,3 g/L	6,6 %
	(197) Analyticon	>>>	1	Negative, < 0,3 g/L	1,6 %
	(201) Iris Diagnostics	>>>	1	Negative, < 0,3 g/L	1,6 %
	(202) Dirui	>>>	19	Negative, < 0,3 g/L	31 %
	(207) 77 Elektronika Kft.	>>>	5	Negative, < 0,3 g/L	8,2 %
	(999) another manufacturer	>>>	1	Negative, < 0,3 g/L	1,6 %
Total per sample: n = 61		Success: 100 %			
B	(5) Dialab	>>>	1	> 2,5 g/L	1,6 %
	(12) Beckman Coulter	>>>	2	0,51 - 2,5 g/L	3,3 %
	(18) Cormay	>>>	4	> 2,5 g/L	6,6 %
	(46) Erba Lachema	>>>	7	0,51 - 2,5 g/L	11 %
		>>>	4	> 2,5 g/L	6,6 %
	(60) Roche	>>>	2	> 2,5 g/L	3,3 %
	(63) Sysmex	>>>	2	0,51 - 2,5 g/L	3,3 %
	(109) Axis-Shield	>>>	2	> 2,5 g/L	3,3 %
	(150) Arkray		1	0,3 - 0,5 g/L	1,6 %
		>>>	1	0,51 - 2,5 g/L	1,6 %

Summary statistics - qualitative results

(Groups: manufacturer of kit)

EQA round: DGP2/22 - Urine Strip Tests Analysis

Deadline: 20.05.2022

Test	Sample	Group	N	Result	N _{rel}	
B	(161) YD Diagnostics	>>>	1	> 2,5 g/L	1,6 %	
	(162) Siemens (Atellica)	>>>	1	0,3 - 0,5 g/L	1,6 %	
	(179) Siemens	>>>	2	0,51 - 2,5 g/L	3,3 %	
	(197) Analyticon	>>>	2	0,51 - 2,5 g/L	3,3 %	
	(201) Iris Diagnostics	>>>	1	> 2,5 g/L	1,6 %	
	(202) Dirui	>>>	1	0,51 - 2,5 g/L	1,6 %	
	(207) 77 Elektronika Kft.	>>>	17	0,51 - 2,5 g/L	28 %	
	(999) another manufacturer	>>>	2	> 2,5 g/L	3,3 %	
	(999) another manufacturer	>>>	3	0,3 - 0,5 g/L	4,9 %	
	(999) another manufacturer	>>>	2	0,51 - 2,5 g/L	3,3 %	
	(999) another manufacturer	>>>	1	> 2,5 g/L	1,6 %	
	Total per sample: n = 61			Success: 92 %		
	Total success: 92 %					

(261) Glucose

A	(5) Dialab	>>>	1	Negative, normal, <= 2,8 mmol/L	1,6 %
	(12) Beckman Coulter	>>>	1	Negative, normal, <= 2,8 mmol/L	1,6 %
	(18) Cormay	>>>	1	2,81 - 11 mmol/L	1,6 %
	(18) Cormay	>>>	1	Negative, normal, <= 2,8 mmol/L	1,6 %
	(18) Cormay	>>>	2	2,81 - 11 mmol/L	3,3 %
	(18) Cormay	>>>	1	11,1 - 28 mmol/L	1,6 %
	(46) Erba Lachema	>>>	8	Negative, normal, <= 2,8 mmol/L	13 %
	(46) Erba Lachema	>>>	3	2,81 - 11 mmol/L	4,9 %
	(60) Roche	>>>	2	Negative, normal, <= 2,8 mmol/L	3,3 %
	(63) Sysmex	>>>	2	Negative, normal, <= 2,8 mmol/L	3,3 %
	(109) Axis-Shield	>>>	2	Negative, normal, <= 2,8 mmol/L	3,3 %
	(150) Arkray	>>>	2	Negative, normal, <= 2,8 mmol/L	3,3 %
	(161) YD Diagnostics	>>>	1	2,81 - 11 mmol/L	1,6 %
	(162) Siemens (Atellica)	>>>	3	Negative, normal, <= 2,8 mmol/L	4,9 %
	(179) Siemens	>>>	4	Negative, normal, <= 2,8 mmol/L	6,6 %
	(197) Analyticon	>>>	1	Negative, normal, <= 2,8 mmol/L	1,6 %
	(201) Iris Diagnostics	>>>	1	Negative, normal, <= 2,8 mmol/L	1,6 %
	(202) Dirui	>>>	19	Negative, normal, <= 2,8 mmol/L	31 %
	(207) 77 Elektronika Kft.	>>>	3	Negative, normal, <= 2,8 mmol/L	4,9 %
	(207) 77 Elektronika Kft.	>>>	2	2,81 - 11 mmol/L	3,3 %
(999) another manufacturer	>>>	1	2,81 - 11 mmol/L	1,6 %	
Total per sample: n = 61			Success: 82 %		
B	(5) Dialab	>>>	1	> 28 mmol/L	1,6 %
	(12) Beckman Coulter	>	1	11,1 - 28 mmol/L	1,6 %
	(18) Cormay	>>>	1	> 28 mmol/L	1,6 %
	(18) Cormay	>>>	1	2,81 - 11 mmol/L	1,6 %
	(18) Cormay	>	1	11,1 - 28 mmol/L	1,6 %
	(18) Cormay	>>>	2	> 28 mmol/L	3,3 %
	(46) Erba Lachema	>>>	1	Negative, normal, <= 2,8 mmol/L	1,6 %
	(46) Erba Lachema	>	7	11,1 - 28 mmol/L	11 %
	(46) Erba Lachema	>>>	3	> 28 mmol/L	4,9 %
	(60) Roche	>>>	2	> 28 mmol/L	3,3 %
	(63) Sysmex	>	1	11,1 - 28 mmol/L	1,6 %
	(109) Axis-Shield	>>>	1	> 28 mmol/L	1,6 %
	(109) Axis-Shield	>>>	2	> 28 mmol/L	3,3 %
	(150) Arkray	>	1	11,1 - 28 mmol/L	1,6 %
	(150) Arkray	>>>	1	> 28 mmol/L	1,6 %
	(161) YD Diagnostics	>>>	1	> 28 mmol/L	1,6 %
	(162) Siemens (Atellica)	>>>	1	2,81 - 11 mmol/L	1,6 %
	(162) Siemens (Atellica)	>>>	2	> 28 mmol/L	3,3 %
	(179) Siemens	>	4	11,1 - 28 mmol/L	6,6 %
	(197) Analyticon	>>>	1	> 28 mmol/L	1,6 %
(201) Iris Diagnostics	>>>	1	> 28 mmol/L	1,6 %	
(202) Dirui	>	6	11,1 - 28 mmol/L	9,8 %	
(202) Dirui	>>>	13	> 28 mmol/L	21 %	
(207) 77 Elektronika Kft.	>	1	11,1 - 28 mmol/L	1,6 %	
(207) 77 Elektronika Kft.	>>>	4	> 28 mmol/L	6,6 %	
(999) another manufacturer	>>>	1	2,81 - 11 mmol/L	1,6 %	
Total per sample: n = 61			Success: 93 %		
Total success: 77 %					

(275) hCG

A	(5) Dialab	>>>	2	Negative	40 %
	(17) Biogema	>>>	1	Negative	20 %

Summary statistics - qualitative results

(Groups: manufacturer of kit)

EQA round: DGP2/22 - Urine Strip Tests Analysis

Deadline: 20.05.2022

Test	Sample	Group	N	Result	N _{rel}
A	(32)	Human	>>>	1 Negative	20 %
	(179)	Siemens	>>>	1 Negative	20 %
		Total per sample: n = 5	Success: 100 %		
B	(5)	Dialab	>>>	2 Positive	40 %
	(17)	Biogema	>>>	1 Positive	20 %
	(32)	Human	>>>	1 Positive	20 %
	(179)	Siemens	>>>	1 Positive	20 %
		Total per sample: n = 5	Success: 100 %		
Total success: 100 %					

(263) Ketone bodies

A	(5)	Dialab	>>>	1 Negative, <= 0,5 mmol/L	1,6 %
	(12)	Beckman Coulter	>>>	2 Negative, <= 0,5 mmol/L	3,3 %
	(18)	Cormay	>>>	3 Negative, <= 0,5 mmol/L	4,9 %
				1 1+, 2+	1,6 %
	(46)	Erba Lachema	>>>	11 Negative, <= 0,5 mmol/L	18 %
	(60)	Roche	>>>	2 Negative, <= 0,5 mmol/L	3,3 %
	(63)	Sysmex	>>>	2 Negative, <= 0,5 mmol/L	3,3 %
	(109)	Axis-Shield	>>>	2 Negative, <= 0,5 mmol/L	3,3 %
	(150)	Arkray	>>>	2 Negative, <= 0,5 mmol/L	3,3 %
	(161)	YD Diagnostics	>>>	1 Negative, <= 0,5 mmol/L	1,6 %
	(162)	Siemens (Atellica)	>>>	3 Negative, <= 0,5 mmol/L	4,9 %
	(179)	Siemens	>>>	4 Negative, <= 0,5 mmol/L	6,6 %
	(197)	Analyticon	>>>	1 Negative, <= 0,5 mmol/L	1,6 %
	(201)	Iris Diagnostics	>>>	1 Negative, <= 0,5 mmol/L	1,6 %
	(202)	Dirui	>>>	19 Negative, <= 0,5 mmol/L	31 %
	(207)	77 Elektronika Kft.	>>>	5 Negative, <= 0,5 mmol/L	8,2 %
(999)	another manufacturer	>>>	1 Negative, <= 0,5 mmol/L	1,6 %	
		Total per sample: n = 61	Success: 98 %		
B	(5)	Dialab	>>>	1 1+, 2+	1,6 %
	(12)	Beckman Coulter	>>>	2 1+, 2+	3,3 %
	(18)	Cormay	>>>	3 1+, 2+	4,9 %
			>	1 3+, 4+	1,6 %
	(46)	Erba Lachema	>>>	8 1+, 2+	13 %
			>	3 3+, 4+	4,9 %
	(60)	Roche	>	2 3+, 4+	3,3 %
	(63)	Sysmex		1 Negative, <= 0,5 mmol/L	1,6 %
			>>>	1 1+, 2+	1,6 %
	(109)	Axis-Shield	>>>	2 1+, 2+	3,3 %
	(150)	Arkray	>	2 3+, 4+	3,3 %
	(161)	YD Diagnostics	>>>	1 1+, 2+	1,6 %
	(162)	Siemens (Atellica)	>>>	3 1+, 2+	4,9 %
	(179)	Siemens	>>>	4 1+, 2+	6,6 %
	(197)	Analyticon	>>>	1 1+, 2+	1,6 %
	(201)	Iris Diagnostics	>>>	1 1+, 2+	1,6 %
(202)	Dirui	>>>	12 1+, 2+	20 %	
		>	7 3+, 4+	11 %	
(207)	77 Elektronika Kft.	>>>	5 1+, 2+	8,2 %	
(999)	another manufacturer	>>>	1 1+, 2+	1,6 %	
		Total per sample: n = 61	Success: 98 %		
Total success: 97 %					

(272) Blood

A	(5)	Dialab	>>>	1 Negative	1,6 %
	(12)	Beckman Coulter	>>>	2 Negative	3,3 %
	(18)	Cormay	>>>	3 Negative	4,9 %
				1 > 50 ery/ μ L	1,6 %
	(46)	Erba Lachema	>>>	11 Negative	18 %
	(60)	Roche	>>>	2 Negative	3,3 %
	(63)	Sysmex	>>>	2 Negative	3,3 %
	(109)	Axis-Shield	>>>	2 Negative	3,3 %
	(150)	Arkray	>>>	2 Negative	3,3 %
	(161)	YD Diagnostics	>>>	1 Negative	1,6 %
	(162)	Siemens (Atellica)	>>>	3 Negative	4,9 %
	(179)	Siemens	>>>	4 Negative	6,6 %
	(197)	Analyticon	>>>	1 Negative	1,6 %
	(201)	Iris Diagnostics	>>>	1 Negative	1,6 %
	(202)	Dirui	>>>	19 Negative	31 %
	(207)	77 Elektronika Kft.	>>>	5 Negative	8,2 %
(999)	another manufacturer	>>>	1 Negative	1,6 %	

Summary statistics - qualitative results

(Groups: manufacturer of kit)

EQA round: DGP2/22 - Urine Strip Tests Analysis

Deadline: 20.05.2022

Test	Sample	Group	N	Result	N _{rel}
Total per sample: n = 61 Success: 98 %					
B	(5)	Dialab	>>>	1 > 50 ery/ μ L	1,6 %
	(12)	Beckman Coulter	>>>	2 > 50 ery/ μ L	3,3 %
	(18)	Cormay		1 Negative	1,6 %
			>>>	3 > 50 ery/ μ L	4,9 %
	(46)	Erba Lachema	>>>	11 > 50 ery/ μ L	18 %
	(60)	Roche	>>>	2 > 50 ery/ μ L	3,3 %
	(63)	Sysmex	>>>	2 > 50 ery/ μ L	3,3 %
	(109)	Axis-Shield	>>>	2 > 50 ery/ μ L	3,3 %
	(150)	Arkray		1 5 - 12 ery/ μ L	1,6 %
			>>>	1 > 50 ery/ μ L	1,6 %
	(161)	YD Diagnostics	>>>	1 > 50 ery/ μ L	1,6 %
	(162)	Siemens (Atellica)	>	1 13 - 50 ery/ μ L	1,6 %
			>>>	2 > 50 ery/ μ L	3,3 %
	(179)	Siemens	>>>	4 > 50 ery/ μ L	6,6 %
	(197)	Analyticon	>>>	1 > 50 ery/ μ L	1,6 %
	(201)	Iris Diagnostics	>>>	1 > 50 ery/ μ L	1,6 %
	(202)	Dirui	>	1 13 - 50 ery/ μ L	1,6 %
			>>>	18 > 50 ery/ μ L	30 %
	(207)	77 Elektronika Kft.	>>>	5 > 50 ery/ μ L	8,2 %
	(999)	another manufacturer	>>>	1 > 50 ery/ μ L	1,6 %
Total per sample: n = 61 Success: 97 %					
Total success: 95 %					

(270) WBC

A	(5)	Dialab	>>>	1 Negative	1,8 %
	(12)	Beckman Coulter	>>>	2 Negative	3,6 %
	(18)	Cormay	>>>	4 Negative	7,1 %
	(46)	Erba Lachema	>>>	8 Negative	14 %
	(60)	Roche	>>>	1 Negative	1,8 %
	(63)	Sysmex	>>>	2 Negative	3,6 %
	(109)	Axis-Shield	>>>	1 Negative	1,8 %
				1 10 - 50 leu/ μ L	1,8 %
	(150)	Arkray	>>>	2 Negative	3,6 %
	(162)	Siemens (Atellica)	>>>	3 Negative	5,4 %
	(179)	Siemens	>>>	4 Negative	7,1 %
	(197)	Analyticon	>>>	1 Negative	1,8 %
	(201)	Iris Diagnostics	>>>	1 Negative	1,8 %
	(202)	Dirui	>>>	19 Negative	34 %
	(207)	77 Elektronika Kft.	>>>	5 Negative	8,9 %
	(999)	another manufacturer	>>>	1 Negative	1,8 %
Total per sample: n = 56 Success: 98 %					
B	(5)	Dialab	>>>	1 > 250 leu/ μ L	1,8 %
	(12)	Beckman Coulter		2 Negative	3,6 %
	(18)	Cormay	>>>	4 51 - 250 leu/ μ L	7,1 %
	(46)	Erba Lachema	>>>	8 > 250 leu/ μ L	14 %
	(60)	Roche	>>>	1 > 250 leu/ μ L	1,8 %
	(63)	Sysmex	>	1 10 - 50 leu/ μ L	1,8 %
			>>>	1 51 - 250 leu/ μ L	1,8 %
	(109)	Axis-Shield	>>>	2 > 250 leu/ μ L	3,6 %
	(150)	Arkray	>>>	1 51 - 250 leu/ μ L	1,8 %
			>>>	1 > 250 leu/ μ L	1,8 %
	(162)	Siemens (Atellica)	>>>	1 51 - 250 leu/ μ L	1,8 %
			>>>	2 > 250 leu/ μ L	3,6 %
	(179)	Siemens	>>>	1 51 - 250 leu/ μ L	1,8 %
			>>>	3 > 250 leu/ μ L	5,4 %
	(197)	Analyticon	>	1 10 - 50 leu/ μ L	1,8 %
	(201)	Iris Diagnostics	>>>	1 > 250 leu/ μ L	1,8 %
	(202)	Dirui	>	1 10 - 50 leu/ μ L	1,8 %
			>>>	17 51 - 250 leu/ μ L	30 %
			>>>	1 > 250 leu/ μ L	1,8 %
	(207)	77 Elektronika Kft.	>>>	1 51 - 250 leu/ μ L	1,8 %
			>>>	4 > 250 leu/ μ L	7,1 %
	(999)	another manufacturer	>>>	1 > 250 leu/ μ L	1,8 %
Total per sample: n = 56 Success: 96 %					
Total success: 95 %					

(274) Nitrite

A	(5)	Dialab	>>>	1 Negative	1,7 %
	(12)	Beckman Coulter	>>>	2 Negative	3,3 %
	(18)	Cormay	>>>	3 Negative	5,0 %

Summary statistics - qualitative results

(Groups: manufacturer of kit)

EQA round: DGP2/22 - Urine Strip Tests Analysis

Deadline: 20.05.2022

Test	Sample	Group	N	Result	N _{rel}
A	(18)	Cormay	1	Positive	1,7 %
	(46)	Erba Lachema	>>> 10	Negative	17 %
	(60)	Roche	>>> 2	Negative	3,3 %
	(63)	Sysmex	>>> 2	Negative	3,3 %
	(109)	Axis-Shield	>>> 2	Negative	3,3 %
	(150)	Arkray	>>> 2	Negative	3,3 %
	(161)	YD Diagnostics	>>> 1	Negative	1,7 %
	(162)	Siemens (Atellica)	>>> 3	Negative	5,0 %
	(179)	Siemens	>>> 4	Negative	6,7 %
	(197)	Analyticon	>>> 1	Negative	1,7 %
	(201)	Iris Diagnostics	>>> 1	Negative	1,7 %
	(202)	Dirui	>>> 19	Negative	32 %
	(207)	77 Elektronika Kft.	>>> 5	Negative	8,3 %
	(999)	another manufacturer	>>> 1	Negative	1,7 %
Total per sample: n = 60			Success: 98 %		
B	(5)	Dialab	>>> 1	Positive	1,7 %
	(12)	Beckman Coulter	>>> 2	Positive	3,3 %
	(18)	Cormay	>>> 4	Positive	6,7 %
	(46)	Erba Lachema	>>> 10	Positive	17 %
	(60)	Roche	>>> 2	Positive	3,3 %
	(63)	Sysmex	>>> 2	Positive	3,3 %
	(109)	Axis-Shield	>>> 2	Positive	3,3 %
	(150)	Arkray	>>> 2	Positive	3,3 %
	(161)	YD Diagnostics	>>> 1	Positive	1,7 %
	(162)	Siemens (Atellica)	>>> 3	Positive	5,0 %
	(179)	Siemens	>>> 4	Positive	6,7 %
	(197)	Analyticon	>>> 1	Positive	1,7 %
	(201)	Iris Diagnostics	>>> 1	Positive	1,7 %
	(202)	Dirui	>>> 19	Positive	32 %
(207)	77 Elektronika Kft.	>>> 5	Positive	8,3 %	
(999)	another manufacturer	>>> 1	Positive	1,7 %	
Total per sample: n = 60			Success: 100 %		
Total success: 98 %					

(273) pH

A	(5)	Dialab	>>> 1	6,5	1,6 %
	(12)	Beckman Coulter	> 2	7,0	3,3 %
	(18)	Cormay	3	5,0	4,9 %
	(46)	Erba Lachema	> 1	6,0	1,6 %
			> 1	5,5	1,6 %
			> 5	6,0	8,2 %
			>>> 4	6,5	6,6 %
	> 1	7,0	1,6 %		
	(60)	Roche	> 2	7,0	3,3 %
	(63)	Sysmex	>>> 2	6,5	3,3 %
	(109)	Axis-Shield	> 1	6,0	1,6 %
	>>> 1	6,5	1,6 %		
	(150)	Arkray	>>> 1	6,5	1,6 %
	> 1	7,0	1,6 %		
	(161)	YD Diagnostics	> 1	7,0	1,6 %
	(162)	Siemens (Atellica)	>>> 1	6,5	1,6 %
	> 2	7,0	3,3 %		
	(179)	Siemens	> 4	7,0	6,6 %
	(197)	Analyticon	> 1	6,0	1,6 %
	(201)	Iris Diagnostics	> 1	6,0	1,6 %
(202)	Dirui	> 8	6,0	13 %	
>>> 11	6,5	18 %			
(207)	77 Elektronika Kft.	>>> 4	6,5	6,6 %	
> 1	7,0	1,6 %			
(999)	another manufacturer	> 1	6,0	1,6 %	
Total per sample: n = 61			Success: 95 %		
B	(5)	Dialab	>>> 1	7,0	1,6 %
	(12)	Beckman Coulter	> 1	6,5	1,6 %
	>>> 1	7,0	1,6 %		
	(18)	Cormay	3	5,0	4,9 %
	> 1	6,5	1,6 %		
	(46)	Erba Lachema	> 3	6,5	4,9 %
	>>> 8	7,0	13 %		
(60)	Roche	>>> 2	7,0	3,3 %	
(63)	Sysmex	>>> 2	7,0	3,3 %	

EQA round: DGP2/22 - Urine Strip Tests Analysis

Deadline: 20.05.2022

Test	Sample	Group	N	Result	N _{rel}
B	(109)	Axis-Shield	>>>	2 7,0	3,3 %
	(150)	Arkray	>>>	2 7,0	3,3 %
	(161)	YD Diagnostics	>	1 6,5	1,6 %
	(162)	Siemens (Atellica)	>>>	3 7,0	4,9 %
	(179)	Siemens	>>>	4 7,0	6,6 %
	(197)	Analyticon	>	1 6,0	1,6 %
	(201)	Iris Diagnostics	>	1 6,5	1,6 %
	(202)	Dirui	>	1 6,0	1,6 %
			>	14 6,5	23 %
			>>>	4 7,0	6,6 %
	(207)	77 Elektronika Kft.	>	4 6,5	6,6 %
			>>>	1 7,0	1,6 %
	(999)	another manufacturer	>>>	1 7,0	1,6 %
Total per sample: n = 61			Success: 95 %		

Total success: 93 %

(266) Urobilinogen

A	(5)	Dialab	>>>	1 Negative, < 17 µmol/L	1,6 %
	(12)	Beckman Coulter	>>>	2 Negative, < 17 µmol/L	3,3 %
	(18)	Cormay	>>>	4 Negative, < 17 µmol/L	6,6 %
	(46)	Erba Lachema	>>>	11 Negative, < 17 µmol/L	18 %
	(60)	Roche	>>>	2 Negative, < 17 µmol/L	3,3 %
	(63)	Sysmex	>>>	2 Negative, < 17 µmol/L	3,3 %
	(109)	Axis-Shield	>>>	2 Negative, < 17 µmol/L	3,3 %
	(150)	Arkray	>>>	2 Negative, < 17 µmol/L	3,3 %
	(161)	YD Diagnostics	>>>	1 Negative, < 17 µmol/L	1,6 %
	(162)	Siemens (Atellica)	>>>	3 Negative, < 17 µmol/L	4,9 %
	(179)	Siemens	>>>	4 Negative, < 17 µmol/L	6,6 %
	(197)	Analyticon	>>>	1 Negative, < 17 µmol/L	1,6 %
	(201)	Iris Diagnostics	>>>	1 Negative, < 17 µmol/L	1,6 %
	(202)	Dirui	>>>	19 Negative, < 17 µmol/L	31 %
	(207)	77 Elektronika Kft.	>>>	5 Negative, < 17 µmol/L	8,2 %
	(999)	another manufacturer	>>>	1 Negative, < 17 µmol/L	1,6 %
Total per sample: n = 61			Success: 100 %		
B	(5)	Dialab	>>>	1 3+, 4+	1,6 %
	(12)	Beckman Coulter	>>>	2 3+, 4+	3,3 %
	(18)	Cormay	>>>	4 3+, 4+	6,6 %
	(46)	Erba Lachema	>>>	11 3+, 4+	18 %
	(60)	Roche	>>>	2 3+, 4+	3,3 %
	(63)	Sysmex	>	2 1+, 2+	3,3 %
	(109)	Axis-Shield	>	1 1+, 2+	1,6 %
			>>>	1 3+, 4+	1,6 %
	(150)	Arkray	>	1 1+, 2+	1,6 %
			>>>	1 3+, 4+	1,6 %
	(161)	YD Diagnostics	>>>	1 3+, 4+	1,6 %
	(162)	Siemens (Atellica)	>	3 1+, 2+	4,9 %
	(179)	Siemens	>	2 1+, 2+	3,3 %
			>>>	2 3+, 4+	3,3 %
	(197)	Analyticon	>>>	1 3+, 4+	1,6 %
	(201)	Iris Diagnostics	>>>	1 3+, 4+	1,6 %
	(202)	Dirui	>	19 1+, 2+	31 %
	(207)	77 Elektronika Kft.	>>>	5 3+, 4+	8,2 %
	(999)	another manufacturer	>>>	1 3+, 4+	1,6 %
Total per sample: n = 61			Success: 100 %		

Total success: 100 %